

AED OPERATOR'S CHECKLIST

AED SERIAL NUMBER: _____

LOCATION OF UNIT: _____

DATE						
Check that Active Status Indicator (ASI) is flashing GREEN (If AED has a 9V battery, replace annually)						
Check overall condition of the unit for: - Damage - Dirt/Contamination - Wear						
Check expiry date on electrodes						
Check expiry date on battery						
Comments:						
Inspected by:						

For replacement pads, batteries, accessories or for assistance with your AED please contact:

AEDSolutions
 #11, 4905 - 102 Ave SE | Calgary, AB | T2C 2X7
 800.260.6362 | www.aedsolutions.ca | customercare@aedsolutions.ca

